

RealTime ready

*Custom qPCR Assays and Panels -
simple online configuration and ordering*

Assay ID **Gene Symbol**

100093	AATF
100095	ABL1
100017	ALFM2
100019	ALFM3
100048	ATF5
100049	AVEN
100052	AZU1
100062	BAG3
100068	BAK1
100079	BCL10
100093	BCL2
100095	BCL2A1
100092	BCL2A1
100104	BCL2A1
100108	BCL2A1
100111	BCL2A1
100113	BCL2A1
100115	BECN1
100120	BFAR
100122	BID
100126	BIK
100131	BIRC2
100135	BIRC3
100139	BIRC5

Assay ID 100095
Gene Symbol BCL2A1
Alias ACC-1, ACC-2, BCL2L5, BFL1, GRS, HSPA1
Description BCL2-related protein A1

Add to Panel

RealTime ready

Function tested qPCR assays for human targets of your choice enable expression profiling of your genes of interest in the context of complex biological pathways or other functionally related groups of genes.

RealTime ready qPCR Assays and Panels 4

RealTime ready Configurator 8

RealTime ready Configuration Guide 15

How to Start 16

Search Assays 17

Single Assay Selection and Ordering 21

Custom Panel Configuration 24

Shopping Cart 31

My Orders 34

Cytoplasm



F-κB
p50/p52
RelA

Ubi

radiation
IκBs

RealTime  ready



RealTime ready

Complex Gene Expression Analysis, Made Simple

■ Maximize efficiency


Select assays for your targets of choice out of gene lists of functionally related gene families or biological pathways. Order them as single assays or configure your own custom panel on a LightCycler® 480 Multiwell Plate.

■ Ensure optimal results

Rely on function tested qPCR assays based on proven Universal ProbeLibrary technology, with detailed bioinformatic background information and novel error-detection features.

■ Get the complete picture

Discover the full functionality of your targets in complex biological pathways, and reveal their influence on processes leading to diseases or disorders.



RealTime ready qPCR Assays and Panels

Gene expression profiling of the human genome on the LightCycler® 480 System has never been so easy

Rapidly and flexibly quantify the expression levels of your targets of interest with ready-to-use RealTime ready qPCR Assays. You can select and order single assays for any human target or configure your custom panel with the selected assays on the free, online RealTime ready Configurator. Alternatively, you can simply rely on our pre-configured RealTime ready Focus Panels with a broad range of selected genes from specific pathways or functional groups.

RealTime ready qPCR Assays

RealTime ready qPCR assays each contain target-specific primers and a Universal ProbeLibrary probe. Each assay is qualified after the design process to function in qPCR on the LightCycler® 480 Instrument using a universal biological sample under standard PCR conditions with the LightCycler® 480 Probes Master. Each assay must comply with the following stringent criteria:

- **Linear dynamic range of at least 3 logs**
- **PCR efficiency of 2.0**
- **R² value of standard curve between 0.99 and 1.00**
- **High amplification specificity, no side products in gel analysis**

RealTime ready assays, like all UPL assays, can be used successfully on any real-time PCR Instrument. For best performance, we recommend using RealTime ready DNA Probes Master.

Triplicate analysis shows almost identical CP values and consistent distances in the dilution series




Figure 3: Dilution series of 50ng to 5 pg of the Human RAN/RANP1 gene, analyzed in triplicates.

No side products detected in agarose gel analysis




Figure 1: 48 RealTime ready Assays of the Human GPCR gene family.

Linearity of amplification over five 10-fold dilution steps with efficiency values of 99.2 to 100 %.




Figure 2: Linearity plot (Cp vs log concentration) over five 10-fold dilution steps for 12 assays from the Human Apoptosis Pathway.

Assay Search and Ordering

Use the Configurator to select and order RealTime ready Catalog Assays (300 reactions of 20 µl each), for human targets of your choice; if no Catalog Assay is available, we will design and test a Designer Assay for you.

Assay Technology

RealTime ready Assays are based on the Universal ProbeLibrary (UPL) technology. UPL probes are short hydrolysis probes, labeled at the 5' end with fluorescein (FAM) and at the 3' end with a dark quencher dye. In order to maintain the specificity and melting temperature that hybridizing qPCR probes require, Locked Nucleic Acids (LNA) are incorporated into the sequence of each UPL probe (Fig. 4). The highly specific RealTime ready Assays can therefore be used under standard PCR conditions on all real-time PCR systems.




Figure 4: UPL probes incorporate LNA nucleosides.

Locked Nucleic Acids are a class of nucleic acid analogs where the ribose ring is “locked” with a methylene bridge connecting the 2'-O atom with the 4'-C atom. The locked ribose conformation enhances base stacking and backbone pre-organization, increasing thermal stability and discriminative power of duplexes. LNAs discriminate single base mismatches under conditions not possible with other nucleic acids.

Assay Design Characteristics

RealTime ready Assays are designed using the established Universal ProbeLibrary ProbeFinder Software, which is based on verified design algorithms. For a particular gene, one unique transcript is chosen for assay design. The selection of this representative transcript is based on annotation and cross references within two of the main public gene/transcript annotation resources, **Ensembl** and **Entrez Gene**. If possible, an intron-spanning assay is designed. This procedure is thought to detect a broad range of splice variants of a certain gene. Almost any gene, protein, transcript, sequence ID, or probe ID may be used as a search term. The annotations are based on the GRCh37 gene build released by the **Genome Reference Consortium**. The “Search Tips” on the Configurator provide a detailed description of individual IDs and search functionality.

.....

Ensembl

(<http://www.ensembl.org>)

.....

Entrez Gene

(<http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>)

.....

Genome Reference Consortium

(<http://www.ncbi.nlm.nih.gov/projects/genome/assembly/grc/index.shtml>)

.....

RealTime ready Custom Panels for the LightCycler® 480 System

Custom panels on LightCycler® 480 Multiwell Plates 96 or 384 are composed of RealTime ready Assays of your choice. A multitude of different layouts (Fig. 5, 6) is available for configuration. Assays are pre-plated and dried-down for a reaction volume of 20 µl on 96-well plates and 10 µl on 384-well plates; simply add sample cDNA and RealTime ready DNA Probes Master for fast and reliable results.

On each panel, up to three wells are reserved for reference gene assays of your choice. In addition, most panel layouts feature a pre-plated error detection concept.

The detailed configuration guide on page 15 shows how to easily populate your panels in a large variety of layouts with the assays of your choice. Alternatively, visit www.realtimeready.roche.com to view a short video demonstration and obtain a wealth of additional information.




Figure 5: RealTime ready Custom Panel 96-16+





Figure 6: RealTime ready Custom Panel 384-64+

Error Detection Feature

1 Reference Genes are mandatory controls on each custom panel and need to be selected during the custom panel configuration process.

2 RT positive control checks for degradation of initial RNA and the quality of the RT step. It consists of three assays, each targeting different positions of the same transcript: 3'-end, center region, and 5'-end.

3 RT minus control detects residual genomic DNA. It consists of two identical assays in two wells: the first well is for the cDNA sample and the second well is for the RNA sample.



RealTime ready Focus Panels – Pre-Configured and Ready-to-Use


RealTime ready Focus Panels are designed for expression profiling of genes from a variety of pathways, functional groups, or gene superfamilies. Ready-to-use, function tested qPCR assays are conveniently pre-plated and dried-down on LightCycler® 480 Multiwell Plates 96 or 384. The only components that need to be added are sample cDNA and LightCycler® 480 Probes Master reaction mix. For each panel, the layout and content can be downloaded from the internet to the Sample Editor of the LightCycler® 480 Software. Each panel contains five reverse-transcription (RT) controls for RNA quality: three RT-positive controls to check the quality of the reverse transcription step, and two RT-minus controls to check for residual genomic DNA. Seven reference gene assays serve as PCR controls and also enable the user to perform relative quantification of the target genes. Select from the following RealTime ready Focus Panels:

- **Human Reference Gene Panel, 96 and 384**
- **Human Cell Cycle Regulation Panel, 96 and 384**
- **Human Apoptosis Panel, 96**
- **Human Apoptosis Panel, 384**
- **Human ABC Transporter Panel, 96 and 384**
- **Human Nuclear Receptor Panel, 96 and 384**



RealTime ready Configurator

Function tested qPCR assays and custom panels are now only a few clicks away. Use this free online configuration and shopping tool to search and order single assays or to configure your own custom panels with assays of your choice.



www.realtimeready.roche.com


Graphical Panel Configuration

The RealTime ready Configurator features an intuitive graphical interface. A large variety of layouts (with or without controls) can be simply filled with assays of your choice. Your changes and edits during the configuration process are automatically stored in your password-protected account at every step; your panel content will be available for download after your purchase is completed.

Search Functions

Assay selection and configuration of custom panels is easy with the comprehensive search functions of the RealTime ready Configurator.


- When you want to order a few single assays or configure layouts with a small number of assays, use the **“Search by Keyword”** option to enter up to five search terms at a time.
- When you want to search for a large number of targets that you have already selected, you can use the **“Batch Upload”** option.
- When you want to study regulation of larger gene groups after knockouts or after cell line treatments that influence certain pathways or signaling cascades, you can use one of two convenient functions, **“Search by Pathway”** or **“Search by Focus Panel”**.



- ➕ Add to Panel
- ✎ Edit Panel
- 🔍 View Panel
- ⚙️ Configure
- 🛒 E-Shop

Search by Keyword


This search function provides you with a variety of options to find assays for up to 5 target genes in one search. You can search by gene name, gene symbol, public accession IDs, biological process, protein function, or a combination of up to 5 keywords.



- Keywords can be single words or short sentences/phrases. Keywords separated by a space are interpreted as exact phrases. Keywords can be refined with the operators AND or OR. For example, 'apoptotic AND inducer' will produce one result for ACIN1 (the 'apoptotic chromatin condensation inducer 1') whereas 'apoptotic inducer' will produce no hits, as the exact phrase is not contained in the description.
- When using both AND and OR in your search, note that 'AND' takes priority over 'OR' in order of operation. For example, 'cullin 1 OR CCND1 AND cell cycle' will be interpreted as
 - Find genes containing 'cullin 1'
 - Find genes containing 'CCND1 AND cell cycle'
- You can limit the search to certain ID or accession number contexts to get faster and more refined results.
- You can use the "wild card" search option to increase your number of hits and maximize the flexibility of your search. For more information, the "Search Tips" in the RealTime ready Configurator provide detailed assistance for using the "Search by Keyword" function.

Search by Focus Panel Lists

With the "Search by Focus Panel" function you can configure your own custom Focus Panel by searching a large number of pre-selected targets of relevant gene superfamilies, such as GPCR, protein kinases, or functionally related groups such as transcription factors or ion channels. This function is also helpful when you just want to modify a RealTime ready Focus Panel with targets of your specific interest.



The tutorial video on www.realtimeready.roche.com walks you through the features and search functions of the RealTime ready Configurator.

The “Search by Pathway” function allows you to select your targets from a variety of interactive maps of cellular or signaling pathways. These maps immediately reveal possible interactions of your selected targets with other genes. Simply click on the genes of interest to add them to your assay search.



Use the “Batch Upload” option to upload up to 384 different assay or accession IDs directly into the assay search process.

A maximum of 384 lines are processed. IDs and accession numbers of various origins may be combined in one file.

You can use accession IDs from the fields listed below for Batch Upload. We recommend using as few different ID types as possible and activating only the appropriate fields for the search in order to get more concise search results (*e.g.*, providing a list of HGNC IDs and applying the filter accordingly will omit any search results where the HGNC ID matches, for example, EntrezGene ID by chance).

| 11

Search Results

In the search results, all genes matching your search term are listed with the gene symbol, aliases, and description, together with the relevant part of your query. If there is already a RealTime ready Assay available for that gene, the transcript accession number that was used for assay design is also provided. For each assay, additional detailed information is provided about the assay and its target, such as the exact location of the assay on the transcript. In case no RealTime ready Catalog Assay is available yet, you can order a Designer Assay for this specific target gene. Assay details will then be provided after the assay has been produced.

Your search by Focus Panel "Focus Panel MD 1" 1-10 of 100 Results Results per page 10

1 | 2 | 3... 10 next >

<input type="checkbox"/> Assay ID	<input type="checkbox"/> Gene Symbol	<input type="checkbox"/> Alias	Description	Transcript	
<input type="checkbox"/> 100003	AATF	CHE-1, CHE1, DED	apoptosis antagonizing transcription factor	ENST00000225402	Details Add to List
<input type="checkbox"/> 100005	ABL1	ABL, bcr/abl, c-ABL, JTK7, p150, v-abl	c-abl oncogene 1, receptor tyrosine kinase	ENST00000372348	Details Add to List
<input type="checkbox"/> 100015	AIFM1	AIF, PDCD8	apoptosis-inducing factor, mitochondrion-associated, 1	ENST00000319908	Details Add to List
<input type="checkbox"/> 100017	AIFM2	AMID, FLJ14497, PRG3	apoptosis-inducing factor, mitochondrion-associated, 2	ENST00000373248	Details Add to List
<input type="checkbox"/> 100019	AIFM3	AIFL, FLJ30473, FLJ45137	apoptosis-inducing factor, mitochondrion-associated, 3	ENST00000399167	Details Add to List
<input type="checkbox"/> 100005	ABL1	ABL, bcr/abl, c-ABL, JTK7, p150, v-abl	c-abl oncogene 1, receptor tyrosine kinase	ENST00000372348	Details Add to List
<input type="checkbox"/> 100015	AIFM1	AIF, PDCD8	apoptosis-inducing factor, mitochondrion-associated, 1	ENST00000319908	Details Add to List
<input type="checkbox"/> 100017	AIFM2	AMID, FLJ14497, PRG3	apoptosis-inducing factor, mitochondrion-associated, 2	ENST00000373248	Details Add to List
<input type="checkbox"/> 100019	AIFM3	AIFL, FLJ30473, FLJ45137	apoptosis-inducing factor, mitochondrion-associated, 3	ENST00000399167	Details Add to List
<input type="checkbox"/> 100005	ABL1	ABL, bcr/abl, c-ABL, JTK7, p150, v-abl	c-abl oncogene 1, receptor tyrosine kinase	ENST00000372348	Details Add to List
<input type="checkbox"/> 100015	AIFM1	AIF, PDCD8	apoptosis-inducing factor, mitochondrion-associated, 1	ENST00000319908	Details Add to List
<input type="checkbox"/> 100017	AIFM2	AMID, FLJ14497, PRG3	apoptosis-inducing factor, mitochondrion-associated, 2	ENST00000373248	Details Add to List
<input type="checkbox"/> 100019	AIFM3	AIFL, FLJ30473, FLJ45137	apoptosis-inducing factor, mitochondrion-associated, 3	ENST00000399167	Details Add to List
<input type="checkbox"/> 100005	ABL1	ABL, bcr/abl, c-ABL, JTK7, p150, v-abl	c-abl oncogene 1, receptor tyrosine kinase	ENST00000372348	Details Add to List
<input type="checkbox"/> 100015	AIFM1	AIF, PDCD8	apoptosis-inducing factor, mitochondrion-associated, 1	ENST00000319908	Details Add to List

Assay Details

For each RealTime ready Assay, the following graphical information is provided:


- Accession ID of the transcript for which the assay was designed
- Postion of the amplicon on this transcript
- Amplicon length
- Intron/exon boundaries for intron-spanning assays

For assay design, the most representative or “main” transcript for an individual gene is identified (“Assay Design Characteristics” on page 5). Therefore, RealTime ready Assays are ideally suited for standard gene expression quantification experiments.

The base-precise location of the assay on other splice variants represented by Ensembl or RefSeq transcripts is shown under the “Transcript Information” on the Assay Details page.

Probe IDs from Affymetrix, Illumina, and NimbleGen gene expression arrays that match any transcripts of a gene for which the assay was designed are listed under “Array Platform IDs”.

Detailed annotation information is also listed for the respective target on the “Assay Details” page.



Gene/Protein Information

Gene and protein IDs associated with any transcript of the gene for which the assay was designed are listed, along with links to the information on the respective databases (if applicable).

▼ Gene / Protein Information	
Ensembl Gene ID	ENSG00000148400
EntrezGene ID	4851
MIM Gene Accession	109730 190198
MIM Morbid Accession	109730
UCSC ID	uc004chz
Unigene ID	Hs.495473
<hr/>	
Ensembl Protein ID	ENSP00000277541 ENSP00000360765

Array Platform IDs

Probe IDs annotated to the respective gene for which the assay was designed are listed, along with links to information on the respective databases.

AFFY_HuGene_1_0_st_v1	8165217
Affy hg u133 plus 2 ID	218902_at 223508_at
Affy hg u133a 2 ID	218902_at
Affy hg u133a ID	218902_at
Affy hg u133b ID	223508_at
Affy hg u95c ID	65751_at
Affy hg u95d ID	86254_r_at
Affy hg u95e ID	76569_at
Affy u133 x3p ID	g8923012_3p_at
Agilent Probe	A_23_P60387 A_23_P60393
Illumina v1 ID	GI_27894367-S

Transcript Information


In concordance with “Gene/Protein Information”, the “Transcript Information” section lists all transcripts associated with the gene for which the assay was designed. If applicable, all listed IDs provide a link to information on the respective databases. Whenever an assay covers an Ensembl or RefSeq transcript, additional precise assay mapping positions are provided. Any potential splice variant specificity of a RealTime ready Assay can thus be deduced.

▼ Transcript Information	
CCDS ID	CCDS43905
EMBL ID	AB209873 AC_000052 AC_000141 AF308602 AK000012 AL354671 AL592301 BC013208 BC039147 BC046127 BC049843 BC063597 CH471090 CN431067 CR457221 DA324222 M73980 NC_000009 NG_007458 NT_024000 NW_001839245 NW_924573
Ensembl Transcript ID	ENST00000361528 ENST00000371700
HAVANA transcript ID (OTTT)	OTTHUMT00000055087
Refseq DNA ID	NM_017617

In addition, the different categories (e.g. Biotype, Codelink, HGNC ID) matching the transcripts are provided as well as the complete GO annotation of the gene detected by the RealTime ready Assay.

RealTime ready Configuration Guide

This tutorial shows you how to search for assays, configure a custom panel, use your custom database, and order your custom products.



How to Start

Logging In

Simply log in with your user ID and password from the Roche Applied Science e-shop **1**, or create a new e-shop account **2**.

Roche Applied Science United States Help Contact Us

RealTime ready Products RealTime ready Special Interest Site

RealTime ready Configurator

1

Registered Users Login

If you are already a registered user of our e-shop, please use your User ID and Password to login. If not, please register below!

> Forgot your Password?

User ID:

Password:

Enter

2

No Account Yet? Register Here!

Registration is quick and free, and allows you to use all functions of the RealTime ready Configurator and the full convenience of online ordering.

Register

Benefits

Registration allows you the convenience of online shopping, the versatility of customizable assay configurations, and the ability to save your work and return to it at any time.

- Enjoy free online storage for all your configurations.
- Retrieve, modify, and update your configurations at any time.
- Your account information is protected; no one but you will see your configurations.

Need Help?

Get in touch with our sales and support teams via phone, fax, or e-mail.

> Contact Us

Visit the RealTime ready Special Interest Site (SIS) or refer to the Configuration Guide:

- > RealTime ready SIS
- > Configuration Guide

© 1996-2009 Roche Diagnostics Corporation Legal Statement Privacy Policy About Us

This website contains information on products which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin.

Logging into your Configurator account is required for online ordering and secure storage of your assay selections and panel configurations.

Selecting the Operating Mode

After logging in, the workflow selection window shows three options: Single Assays **3**, Custom Panels **4**, and My Configurations **5**.

Roche Applied Science United States

RealTime ready Configurator

3

Single Assays

Go directly to the assay search tool.

Enter

4

Custom Panels

Use the online design tool to build and edit your own panel.

Enter

5

My Configurations

View your saved configurations, make changes to them, and order online, via fax, or e-mail.

Enter


Last Name |

© 1996-2009 Roche Diagnostics Corporation Legal Statement Privacy Policy About Us


This website contains information on products which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin.

Search Assays

Select a Search Function



Selecting “Single Assays” 1 leads you directly to the “Search Assays” page. Choose your preferred search function by clicking tabs 2 - 5.



Search by Keyword

This search function allows you to enter up to five search terms in any combination from any or all of the following five categories: gene symbol, gene name, public accession ID number, biological process, or molecular function.



The screenshot shows the 'Search Assays' interface. At the top, there's a header with instructions: 'Use the different search functions to select Catalog Assays or have us design and test a Designer Assay of your choice. Assays are supplied with 300 x 20ul reactions in solution. For best results, use the RealTime ready Probes Master in your PCR reaction.' Below this are five tabs: 'Search by Keyword', 'Search by Focus Panel', 'Batch Upload', 'Search by Pathway', and 'Reference Genes'. The 'Search by Keyword' tab is active. It contains a search box with the text 'cullin 1 OR CCND1 AND cell cycle'. To the left of the search box are two checkboxes: 'Wildcard search disabled' (checked) and 'Search all Assays' (unchecked). Below these is a 'Limit search to:' dropdown menu. To the right of the search box is a 'Reset' button and a 'Search' button. On the far right, there's a 'My Assays' section with 'Catalog Assays: 0' and 'Designer Assays: 0'. Below that is a 'Search Tips' box with the text 'Could not find what you were looking for? Look at our suggestions to optimize your search.' and a link '> Search Tips'. Numbered callouts are placed over the interface: 1 points to the 'Wildcard search disabled' checkbox, 2 points to the 'Search all Assays' checkbox, 3 points to the 'Search' button, and 4 points to the 'Search Tips' box.

When you want to limit your search to assays for a gene containing all the keywords, type “AND”, in full capitals, between each search term (e.g., “CCND1 AND cell cycle”).

When you want to perform a broader search for all genes that contain any of the keywords, type “OR”, again in full capitals, between the search terms (e.g., “CCND1 OR cell cycle”).


Use the asterisk character as a ‘wildcard’, or placeholder character, to produce a more general search. When using e.g., CCN*, you will find all gene names beginning with the keyword CCN. Note that wildcard searching is disabled by default.

To activate wildcard searching, uncheck this box **1**.

To speed up your search, you can limit the search to different search fields **2**.

Start your search by clicking the “Search” tab **3**. If you have difficulties, use the “Search Tips” box **4** to get helpful information about the search fields.

Search by Focus Panel



Search by Keyword **1 Search by Focus Panel** Batch Upload Search by Pathway


Choose your assays from pre-generated panels of specific pathways, gene/protein families or functional groups. Search for a Focus Panel by the first letter of the name.

1-9 | A | B | C | D | E | F | G | H | I | J | K | L | M | N | O | **P** | Q | R | S | T | U | V | W | X | Y | Z | # | all

Panel (Description)	Number of Assays included	
PDGF signaling	56	Show Assays
Protease	470	Show Assays
Protease Inhibitors	18	Show Assays
SERPINS		
Protein Kinases	474	Show Assays
List of all proteins with protein kinase activity		
Protein Kinases - AGC group	50	Show Assays
Protein Kinases - Atypical group	31	Show Assays

The “Search by Focus Panel” function **1** allows you to search a large number of pre-selected targets of relevant gene superfamilies, such as GPCR, protein kinases, or functionally related groups, such as transcription factors and ion channels.

Batch Upload



Search by Keyword Search by Focus Panel **2 Batch Upload** Search by Pathway Reference Genes

Use this window to upload a text file, containing up to 384 RealTime ready assay IDs or accession IDs from the fields defined below. You can speed up the search by activating the appropriate field below.

Your upload file should be saved as a text file (.txt) without extra formatting. For best results, separate each search term with a carriage return; do not combine terms on a single line or use punctuation to separate terms.

[+ Upload File](#)

ID Fields

☐ Assay ID

☐ Affymetrix

☐ IPI

☐ Agilent

☐ MIM

☐ CCDS

☐ miRBase

☐ EMBL / Genbank / DDBJ

☐ NimbleGen

☐ Ensembl

☐ PDB

☐ Entrez Gene

☐ Refseq

☐ HAVANA

☐ Rfam

☐ HGNC

☐ UCSC

☐ Illumina

☐ Unigene

☐ IMGT

☐ UniProt / Swissprot / TrEMBL

[Reset](#)

The “Batch Upload” function **2** allows you to upload a text file (.txt) that specifies your targets in the form of accession ID's from defined fields, or RealTime ready Assay IDs. Separate each search term with a carriage return. Do not combine terms on a single line, or use punctuation to separate terms.

Search by Pathway


Search by KeywordSearch by Focus PanelBatch Upload1Search by Pathway

Use our convenient pathway maps to search for assays. Simply click on your targets of choice to add them to the assay search function. Search for a pathway by the first letter of the name.

1-9 | A | B | C | D | E | F | G | H | I | J | K | L | M | N | O | P | Q | R | S | T | U | V | W | X | Y | Z | # | all

Pathway	Description	Pathway Map
NFκB Signaling	Members of the NFκB family of transcription factors are important regulators of cell proliferation and cell survival e.g. by regulation of anti-apoptotic genes like TARF1/2. In addition NFκB has a pivotal role in regulating innate and adaptive immune responses.	Pathway Map
Notch Signaling	The evolutionary conserved Delta/Notch signaling is utilized in many cell fate decisions during developmental processes. The Notch protein life cycle includes 3 distinct proteolytic cleavages that finally lead to the Notch IntraCellular Domain (NICD) which carries out its regulatory functions as a co-activator of CLS.	Pathway Map
Sonic Hedgehog	Hedgehog development in pattern maintain organism	
TNF	TNF receptor signaling	
wnt/frizzled signaling	Signaling of the Wnt receptors morphogenesis development	

The “Search by Pathway” function 1 allows you to select your targets from a variety of interactive maps of cellular or signaling pathways. Simply click on a target 2 to add it to your search.



Single Assay Selection and Ordering

Search and Choose Assays

Depending on the type of search you have selected, you will get a list of assays matched to your search terms. If no matching Catalog Assay is available, the result is marked with "tbd", which means Roche will design and test an appropriate RealTime ready Assay for you.

Catalog Assays are accompanied by detailed information. For any displayed assay, the "Details" **1** or "Assay ID" **2** links provide you with an overview of the assay details, including primer-probe location on the transcript, amplicon and transcript length, and intron/exon spanning positions. For a detailed description, please see page 13.

The screenshot shows a search results page for the keyword "cullin 1". The results table lists assays with columns for Assay ID, Gene Symbol, Alias, Description, Transcript, and Search Term. Numbered callouts indicate: **2** for the Assay ID link, **3** for the sort arrows, **4** for the Transcript link, and **1** for the Details link. Below the table, a "RealTime ready Assay Details" panel is shown for Assay ID 100309, displaying a genomic map and assay parameters.

Assay ID	Gene Symbol	Alias	Description	Transcript	Search Term
101454	CUL1		cullin 1	ENST00000325222	cullin 1
100309	CUL1		cullin 1	ENST00000325222	cullin 1
tbd	DCUN1D1	DCUN1L1, RP42, SCCRO, SCRO, Tes3	DCN1, defective in cullin neddylation 1, domain containing 1 (S. cerevisiae)	cullin 1	
tbd	DCUN1D2	C13orf17, FLJ10704, FLJ20092	DCN1, defective in cullin neddylation 1, domain containing 2 (S. cerevisiae)	cullin 1	
tbd	DCUN1D3	44M2, DKFZp686O0290, FLJ41725, MGC48972	DCN1, defective in cullin neddylation 1, domain containing 3 (S. cerevisiae)	cullin 1	
tbd	DCUN1D4	FLJ42355, KIAA0276	DCN1, defective in cullin neddylation 1, domain containing 4 (S. cerevisiae)	cullin 1	
tbd	DCUN1D5	FLJ32431, MGC2714	DCN1, defective in cullin neddylation 1, domain containing 5 (S. cerevisiae)	cullin 1	

RealTime ready Assay Details

Assay ID: 100309, Gene Symbol: CUL1, Alias: cullin 1, Accession ID: ENST00000325222

Amplimer Length: 79 bp, Intron Spanning / Position: 1097-1175, Transcript Length: 3208 bp, Description: cullin 1, HGNC (curated): CUL1

Click on the small arrows **3** to sort the content of these tabs. Clicking on the Transcript **4** leads directly to the Ensembl database for more detailed information.

The screenshot shows the Ensembl genome browser interface for the CUL1 gene. The "Transcript summary" tab is selected, displaying a genomic map of the CUL1 gene with exons and introns. The transcript CUL1-201 is highlighted. The interface includes a sidebar with navigation options, a main content area with transcript details, and a bottom section with statistics and alternative transcripts.

Transcript summary

Transcript: CUL1-201 (ENST00000325222)

Location: Chromosome 7: 148,395,933-148,398,201 Forward strand

Gene: CUL1

This transcript is a product of gene ENSG00000333330 - There are 3 transcripts in this gene: [hide transcripts](#)

Name	Transcript ID	Protein ID	Description
CUL1-201	ENST00000325222	ENSP00000325222	protein_coding
CUL1-202	ENST00000325222	ENSP00000325222	protein_coding
CUL1-203	ENST00000325222	ENSP00000325222	protein_coding

Transcript and Gene level displays

In Ensembl a gene is made up of one or more transcripts. Views in Ensembl are separated into Gene based views and Transcript based views according to which level the information is more appropriately associated. This view is a transcript level view. To flip between the two sets of views you can click on the Gene and Transcript tabs in the menu bar at the top of the page.

Transcript summary

Statistics: Exons: 22, Transcript length: 3,208 bp, Translation length: 776 residues, Type: Known protein coding

Prediction Method Transcript: where the Ensembl genebuild transcript and the vega manual annotation have the same sequence, for every base pair. See [ENST00000325222](#)

Alternative transcripts This Ensembl/RefSeq merge transcript entry corresponds to the following database identifiers: [ENST00000325222](#) [ENST00000325222](#) [ENST00000325222](#)

Select the desired assay with the “Add to List” button on the appropriate line **5**. To select several assays at once, check the boxes next to the desired assays, or select all assays in the displayed table by clicking the box at the top **6**. Then use the “Add to List” **7** button at the top of the page to transfer all the selected assays to your shopping list.


To edit and review your current selections, click the “Edit List” button **8**.

The screenshot shows a search results page for the keyword "cullin 1" in "All". It displays 7 results. Numbered callouts highlight the following features:

- 6**: A checkbox in the top left corner of the table, used to select all displayed assays.
- 5**: The "Add to List" button located at the end of a specific assay row.
- 7**: The "Add to List" button located at the top right of the page, above the table.
- 8**: The "Edit List" button located at the top right of the page, next to the "Add to List" button.

Assay ID	Gene Symbol	Alias	Description	Transcript	Search Term	Actions
<input type="checkbox"/> 101454	CUL1		cullin 1	ENST00000325222	cullin 1	Details Add to List
<input type="checkbox"/> 100309	CUL1		cullin 1	ENST00000325222	cullin 1	Details Add to List
<input type="checkbox"/> tbd	DCUN1D1	DCUN1L1, RP42, SCCRO, SCRO, Tes3	DCN1, defective in cullin neddylation 1, domain containing 1 (S. cerevisiae)		cullin 1	Details Add to List
<input type="checkbox"/> tbd	DCUN1D2	C13orf17, FLJ10704, FLJ20992	DCN1, defective in cullin neddylation 1, domain containing 2 (S. cerevisiae)		cullin 1	Details Add to List
<input type="checkbox"/> tbd	DCUN1D3	44M2, DHTZp666O0290, FLJ41725, MGC48972	DCN1, defective in cullin neddylation 1, domain containing 3 (S. cerevisiae)		cullin 1	Details Add to List
<input type="checkbox"/> tbd	DCUN1D4	FLJ42355, KIAA0276	DCN1, defective in cullin neddylation 1, domain containing 4 (S. cerevisiae)		cullin 1	Details Add to List
<input type="checkbox"/> tbd	DCUN1D5	FLJ32431, MGC2714	DCN1, defective in cullin neddylation 1, domain containing 5 (S. cerevisiae)		cullin 1	Details Add to List


Edit List



The “Edit List” window shows your current assay selections. Moving the mouse cursor over an assay will pop up a box listing brief details for the assay **1**. Use these buttons to edit, refine or print your search results **2**. Note that assays deleted from the “Edit List” page are still available on your “My Configurations” page, in case you want to use them in future. Once you have decided which assays to purchase, move on to the “Review” page **3** to see the list price of your selections.

Review

The “Review” window displays a summary of your selections, along with the list price.




To make changes to your selections, click “Edit List” **1** to return to your list of assays. If you are satisfied with your selections, click “Add to Shopping Cart” **2** to place your assay in the Configurator’s shopping cart in preparation for ordering.

Custom Panel Configuration

The Custom Panels mode **1** allows you to configure RealTime ready Custom Panels with assays of your choice on LightCycler® 480 Multiwell Plates 96 or 384.

The Custom Panel Configurator guides you through this six-step process.



Select Layout and Assign Name

Roche Applied Science | United States | Log Off | Help | Contact Us

RealTime ready Products | RealTime ready Special Interest Site | US First Name US Last Name | empty

RealTime ready Custom Panels

Select Layout | Configuration | Search and Choose Assays | Choose Reference Genes | Review | Order | My Configurations | My Orders

Choose a plate format and a layout type
Choose a plate format and layout, and assign it a unique name which you can use to retrieve, edit, or reorder your custom configuration. You can also add a personalized description.

Name: *
Description:


How to...
Select your plate format and a layout using the "View Layout" function, then click "Configure" to start adding assays to your panel. Please note, that on each panel, you may add up to 8 different designer assays. When you are finished choosing assays and reference genes, you can finalize and order your configuration. Your configurations are safely stored under "My Configurations" at any time and can be edited. Finalized panels cannot be changed, but can be re-ordered and assigned as a...

Layout	Format	Targets	Replicates	Control Wells	Pack Size
LightCycler® 480 Multiwell Plate 96, white					
<input type="checkbox"/> RealTime ready Custom Panel 8*12	96	6	12	no	4
<input type="checkbox"/> RealTime ready Custom Panel 16*6+	96	8	6	yes	
<input type="checkbox"/> RealTime ready Custom Panel 24*4	96	21	4	no	
<input type="checkbox"/> RealTime ready Custom Panel 48*2+	96	40	2	yes	10
<input type="checkbox"/> RealTime ready Custom Panel 96*1	96	93	1		
LightCycler® 480 Multiwell Plate 384, white					
<input type="checkbox"/> RealTime ready Custom Panel 16*24	384	13	24	no	6
<input type="checkbox"/> RealTime ready Custom Panel 32*12+	384	24	12	yes	10

Layouts without control wells contain no controls, but a certain number of reference gene assays must be selected. These reference gene assays can be used as functional controls and as normalization standards in relative quantification analysis of the results.

Visit the RealTime ready Special Interest Site (SIS) or refer to the Configuration Guide:
 > RealTime ready SIS
 > Configuration Guide

Several different layouts are provided for LightCycler® 480 Multiwell Plates. The table provides basic information about each layout. Move the mouse cursor over an entry in the "Control Wells" column **2** to pop up a box containing information about that layout's control concept.



Click on the “View Layout” link to view graphical representations of the panel layouts, and choose between panels with **1** or without **2** controls. Note that all panel layouts require two or three reference genes, which you need to select during the configuration process, described later in this tutorial.

RealTime ready Custom Panels

Select Layout | Configuration | Search and Choose Assays | Choose Reference Genes | Review | Order

My Configurations | My Orders

Choose a plate format and a layout type

Choose a plate format and layout, and assign it a unique name which you can use to retrieve, edit, or reorder your custom configuration. You can also add a personalized description.

3 Name: * Description:

LightCycler® 480 Multiwell Plate 96, white

Layout	Format	Targets	Replicates	Control Wells	Pack Size	
<input type="checkbox"/> RealTime ready Custom Panel 8*12	96	8	12	no	4	View Layout
<input checked="" type="checkbox"/> RealTime ready Custom Panel 16*8+	96	8	6	yes	10	View Layout
<input type="checkbox"/> RealTime ready Custom Panel 24*4	96	21	4	no	10	View Layout
<input type="checkbox"/> RealTime ready Custom Panel 48*2+	96	40	2	yes	10	View Layout
<input type="checkbox"/> RealTime ready Custom Panel 96*1	96	93	1	no	10	View Layout

LightCycler® 480 Multiwell Plate 384, white

5 [Configure](#)


4 Need Help? Get in touch with our sales and support teams via phone, fax, or e-mail. [Contact Us](#)

Visit the RealTime ready Special Interest Site (SIS) or refer to the Configuration Guide: [RealTime ready SIS Configuration Guide](#)


Once you have decided on format and layout, you must provide a unique name for your panel **3**. You may also add a short description for your own records. You can always refer to the “Help” box **4** for more detailed information, technical background, or customer support.

Click “Configure” **5** to start the configuration process.


Configuration



This page gives you a short overview of the configuration workflow **1**, introducing the most important action buttons and icons. The four search functions are also briefly explained. Choosing a search function will take you to the assay search and selection process, as described on pages 17 - 20 **2**.




The “My Panel” window **3** will accompany you throughout the configuration process. It provides you with real-time information about the panel you are currently configuring, including the number of assigned and empty wells, and the number of assay and reference genes.



After finding the assays you need, you can add them to your custom panel. You can select a single assay **4**, or several assays at once by checking multiple boxes **5**, or you can use this top checkbox to mark all assays **6** in the table. Once your assays are selected, click “Add to Panel” **7**.

You can choose up to eight Designer Assays per panel. If you require more than eight Designer Assays or a larger number of panels, please contact your local Roche sales representative.



Use the “View Panel” button **8** to obtain a preview of your custom panel. In the preview window, you can move the mouse cursor over the assays to obtain brief information about each assay. Click the assay to view full information about the assay.

Click ‘Edit Panel’ **9** to work with your current custom panel.

Edit Panel

The “Edit Panel” page shows your current panel content as a list. Move the mouse cursor over the assays to pop up a window with information about the search term and the assay ID **1**.

RealTime ready Custom Panels

Select Layout **Configuration** Search and Choose Assays Choose Reference Genes Review Order My Configurations My Orders

Edit Your Configuration

Below is a list of all the assays you have selected for your panel. Use the mouse over function to view a short description of the assay. You can delete or replace assays or start a new search.

Assay ID	Gene Symbol	Assay ID	Gene Symbol
101464	CUL1	101460	CUL3
100309		100317	CUL4A
101456	CUL1	100320	CUL5
102864	Alias	119315	CUL7

My Panel

Name: Test123456
Layout: RealTime ready Custom Panel 16*6+
Assigned Wells: 8
Empty Wells: 0
Designer Assays: 0
Unique Assays: 8
Reference Genes: 0/3

View Panel
Delete
New Search
Search Results
Choose Reference Genes
Print Panel as List
Review

Before you move on to the “Review” page, select your preferred reference genes **2** for the pre-determined locations on your panel. As you select reference genes, the “My Panel” information box will update accordingly. You can fill additional wells of your panel with reference gene assays of your choice or designate any custom assay as such.

RealTime ready Custom Panels

Select Layout Configuration Search and Choose Assays **Choose Reference Genes** Review Order My Configurations My Orders

Choose Reference Genes

Use this page to select the reference genes for the pre-defined wells on your chosen layout. You can add additional reference genes to your panel by using the “Search by Keyword” function.

Add to Panel View Panel Edit Panel Review

Assay ID	Gene Symbol	Alias	Description	Transcript
102108	ALAS1	ALAS, ALAS3, ALASH, MIG4	aminolevulinate, delta-, synthase 1	ENST00000394965
122801	Ens_G_TSSA0824-00001_symbol	HGNC_automatic_gene_name_TSSA0824-00005, HGNC_curated_gene_name_TSSA0824-00006, HGNC_symbol_TSSA0824-00004	Ens_G_TSSA0824-00001_name	Ens_T_TSSA0824-00000
122802	Ens_G_TSSA0824-00113_symbol	HGNC_automatic_gene_name_TSSA0824-00117, HGNC_curated_gene_name_TSSA0824-00118, HGNC_symbol_TSSA0824-00116	Ens_G_TSSA0824-00113_name	Ens_T_TSSA0824-00172
122803	Ens_G_TSSA0824-00225_symbol	HGNC_automatic_gene_name_TSSA0824-00229, HGNC_curated_gene_name_TSSA0824-00230, HGNC_symbol_TSSA0824-00228	Ens_G_TSSA0824-00225_name	Ens_T_TSSA0824-00284

My Panel

Name: Test09
Layout: RealTime ready Custom Panel 24*4
Assigned Wells: 0
Empty Wells: 21
Designer Assays: 0
Unique Assays: 0
Reference Genes: 0/3

When you are ready, move on to the “Review” page **3**.

Review

The “Review” window summarizes all relevant information about your panel and also provides the price of the panel with all the assays currently selected.

RealTime ready Custom Panels

Select Layout | Configuration | Search and Choose Assays | Choose Reference Genes | **Review** | Order | My Configurations | My Orders

Configuration Review

Please review your configuration. A unique configuration number is assigned when you click "Finalize" or "Add to Shopping Cart".
After we have processed your order, we will provide the panel content in a text file on the "My Orders" page, which can be imported to the sample editor of the LightCycler® 480 Software.

Layout:	RealTime ready Custom Panel 16*6	Name:	Test5
Assigned wells:	3	Description:	for example
Empty wells:	5	Selected Reference Genes:	
Replicates:	6	Price:	452.00 USD
Catalog Assays:	2		
Designer Assays:	1		
Unique Assays:	3		
Control Wells:	yes		
Reference Genes:	0 / 3		
Pack Size:	10 Plates		

[View Panel](#) [Edit Panel](#) [Finalize](#) [Add to Shopping Cart](#)

Need Help?
Get in touch with our sales and support teams via phone, fax, or e-mail.
[Contact Us](#)
Visit the RealTime ready Special Interest Site (SIS) or refer to the Configuration Guide:
[RealTime ready SIS](#)
[Configuration Guide](#)

If you have selected too many or not enough assays or reference genes, helpful reminders will pop up, directing you back to the “Edit Panel” **1** page to make corrections.

To make your purchase immediately, transfer your completed custom panel to the Configurator’s shopping cart **3**.

To place your order later, click “Finalize” **2** to save your panel, with a unique configuration number, in your account’s “My Configurations” area.

My Configurations

Your account's personalized storage and work space, called “My Configurations”, is accessible from the operating mode selection page after logging in, as well as from most other Configurator windows.

The screenshot shows the 'RealTime ready Custom Panels' interface. At the top, there are tabs for 'Select Layout', 'Configuration', 'Search and Choose Assays', 'Choose Reference Genes', 'Review', and 'Order'. A 'My Configurations' button is circled in blue. Below the tabs, there are 'Delete' and 'Add to Shopping Cart' buttons. A text block explains that configurations are stored here and can be viewed, edited, or re-ordered. Below this, there are two sections: 'My Custom Panels' and 'My Single Assays'.

My Custom Panels

Name	Description	Configuration No.	Status	Last Modified	Actions
Test5 RealTime ready Custom Panel 1676	for example		not finalized	2009-09-03	<ul style="list-style-type: none"> Edit Configuration View Panel Export Configuration Delete
Test123456 RealTime ready Custom Panel 364*1			not finalized	2009-08-31	<ul style="list-style-type: none"> Edit Configuration View Panel Export Configuration Delete
Test_04 RealTime ready Custom Panel 6*12		000000727	finalized	2009-09-03	<ul style="list-style-type: none"> Make a Copy View Panel Export Configuration Add to Shopping Cart
Test001 RealTime ready Custom Panel 364*1	Description	000000019	finalized	2009-08-19	<ul style="list-style-type: none"> Make a Copy View Panel Export Configuration Add to Shopping Cart

My Single Assays

Assay ID	Gene Symbol (HGNC)	Alias	Configuration No.	Status	Actions
Designer Assay for AC026283.2-9 Assay ID: 184	AC026283.2-9			not finalized	<ul style="list-style-type: none"> Delete Add to List
Catalog Assay for BCL2 Assay ID: 100083	BCL2	Bcl-2		not finalized	<ul style="list-style-type: none"> Assay Details Delete Add to List
Catalog Assay for CCNH Assay ID: 101394	CCNH	CAK, p34, p37	000000602	finalized	<ul style="list-style-type: none"> Assay Details Add to Shopping Cart

A callout menu is shown on the right, listing the actions available for each configuration: Edit Configuration, View Panel, Export Configuration, and Delete. The 'Delete' button is highlighted in red.

By default, all your configurations and assays are automatically stored as “not finalized”, meaning that you can freely edit and modify them, as well as export them as text files. You can log out of the Configurator, return at a later time, and access all your work in your “My Configurations” page. To continue working on a configuration, click on “Edit Panel” to go to the “Configuration” page, from where you can access all other functions. Click on “Add to List” next to an un-finalized assay to go directly to the “Edit List” page.

Finalized panels and assays have a configuration number assigned to them and cannot be edited or deleted. Finalized panels can, however, be copied and renamed; the duplicate panel can then be edited under the new name. Your original configuration will still be saved under its original name and configuration number, and can be directly placed into the Configurator’s shopping cart.

Shopping Cart

The RealTime ready Configurator's shopping cart shows all custom panels or single assays that you have selected for purchase.

RealTime ready Shopping Cart

Select Layout | Configuration | Search and Choose Assays | Choose Reference Genes | Review | Order | My Configurations | My Orders

Shopping Cart

Please review your purchase selections and quantities. You can continue shopping if you like, or click "E-Mail" when you want to order by e-mail. Use the "Print" function when you want to order by fax. Click "E-Shop" to complete your purchase online via the e-shop. An estimated shipping time (EST) and the country specific catalog price are indicated in the list below. Please note that shipping times will vary depending on the product; you may receive some items before others.

After we have received your order, you can track the actual order status under "My Orders", and download the final panel content as a text file for import into the sample editor of the LightCycler® 480 Software, after your order has been processed.

If you have any questions, please contact us.

Additional terms may apply to custom products and such terms will appear in the Terms and Conditions, quotation or other documentation.

Need Help?
Get in touch with our sales and support teams via phone, fax, or e-mail.
> Contact Us
Visit the RealTime ready Special Interest Site (SIS) or refer to the Configuration Guide:
> RealTime ready SIS
> Configuration Guide

Custom Panels

Configuration No.	Catalog No.	Quantity	Pack Size	Price	Subtotal
000000875	03300228001	1	10 Plates	452.00 USD	452.00 USD

RealTime ready Custom Panel 16*
EST: 18 days

Single Assays

Configuration No.	Catalog No.	Quantity	Price	Subtotal
000000899	04738314001	1	99.00 USD	99.00 USD
000000596	04738314001	1	99.00 USD	99.00 USD

Catalog Assay for CCND1
Assay ID: 100844
300 reactions/20µl
EST: 10 days

Catalog Assay for CCND2
Assay ID: 101384
300 reactions/20µl
EST: 10 days

Products Chosen

Catalog No.	Quantity	Price	Subtotal
04738314001	1	450.00 USD	450.00 USD

Total Price: 452.00 USD

Update Price | Print | E-Mail | E-Shop

Assay Details
The panel consists of four identical sets of pre-tested, ready-to-use qPCR reference gene assays, supplied in a LightCycler® 480 Multiwell Plate (96-well). Each assay includes appropriate primers and a short hydrolysis probe that

You will see that a configuration number **1** has been assigned to each single assay or custom panel. Be sure to use this unique identifying number to ensure accuracy when ordering, since Roche identifies your panels and assays only via this configuration number. An estimated shipping time is provided for each custom product.

You can use the "Assay Details" **3** and "View Panel" **2** options to double-check the contents and details of your panel and assays.

Changes can be made to the number of assays or panels. After making any changes, please remember to click "Update Price" **4** to obtain the new totals for the assays and panels.

Additional Roche Applied Science products, such as the RealTime ready DNA Probes Master, can also be added to the shopping list, to ensure successful performance of your custom panels or assays.

When you are satisfied with the contents of your shopping cart, there are a number of different options to complete your order.


- Click on “E-Mail” **6** and complete the provided order form.
Clicking again on “E-Mail” will electronically send this order to your local sales organization of Roche Applied Science.
- Click “Print” **5** to produce a hard copy of the order form.
Fill out the order form that appears in the window, and then click the “Continue” button, followed by the “Print” button. You can then fax this hard copy to Roche Applied Science to complete your order.
- Click “E-Shop” to place your order online at the Roche Applied Science E-Shop **7**. Using this option will transfer your Configurator’s shopping cart automatically to the shopping cart in the E-Shop.

You will receive a prompt that your Configurator’s shopping cart has been transferred to the Roche Applied Science E-Shop, and a link to the Roche Applied Science home page will appear.

The contents of your shopping cart was transferred to the e-shop. To finalize and submit your order please log-in at roche-applied-science.com and go to the shopping cart.

Click the new link to go to the Roche Applied Science E-Shop.

Roche Applied Science E-shop



When you arrive on the Roche Applied Science home page, please enter your E-Shop User ID and Password **1** and then go to your shopping cart **2**.

Roche Applied Science United States Log Off Quick Order Shopping Cart Help Contact Us

Home Products Special Interest Sites Support & Resources News Search ...

Home > Shopping Cart

Shopping Cart of US First Name US Last Name

3 **1** Successfully updated shopping cart

You have the following items in your Shopping Cart

- If necessary, please adjust the order quantity (change the number)
- Remove items from this order (by placing a check in the Remove box)
- Hit the 'Update the Shopping Cart' button.

Catalog #	Product Name / Pack Size	Configuration #	Price	Qty	Subtotal	Options
03315843001	TxEsXt RTR Single Assay (Designer) 1 test pack	000007816	\$258.00	1	\$258.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
03115856001	TxEsXt RTR Custom Panel 128x3 1 test pack	000004941	\$42.00	1	\$42.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
04738314001	TxEsXt RTR Single Assay (Catalogue) 1 test pack	000007898	\$59.00	2	\$118.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
04738314001	TxEsXt RTR Single Assay (Catalogue) 1 test pack	000007904	\$59.00	1	\$59.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
03300226001	TxEsXt RTR Custom Panel 16x5+ 1 test pack	000007913	\$452.00	1	\$452.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
04738314001	TxEsXt RTR Single Assay (Catalogue) 1 test pack	000007898	\$59.00	2	\$118.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
04738314001	TxEsXt RTR Single Assay (Catalogue) 1 test pack	000007995	\$59.00	2	\$118.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
04738314001	TxEsXt RTR Single Assay (Catalogue) 1 test pack	000007977	\$59.00	1	\$59.00	<input type="checkbox"/> Remove <input type="checkbox"/> Add to favorites
					Total Price	\$1,074.00

Enter a promotion / Quote

4 Update Shopping Cart **5** Continue With Checkout

This green message window **3** should verify that your shopping cart was successfully updated with the contents of your Configurator's shopping cart contents.

You can also delete assays and panels within the E-Shop shopping cart. Be certain to click on "Update Shopping Cart" **4** after making any changes.

After updating and verifying the contents of the E-Shop shopping cart, click "Continue With Checkout" **5** to send your order to Roche.

My Orders

When your custom orders are active and arrive at Roche, you will find them listed in the RealTime ready Configurator under “My Orders” **1**.

RealTime ready Configurator

My Configurator **1** My Orders

Add to Shopping Cart

After we have received your RealTime ready order, you can track the status on this page. If you have ordered a custom panel, we will provide the panel content in a text file, which can be imported to the sample editor of the LightCycler® 480 Software.

▼ My Ordered Custom Panels

Name	Description	Configuration No.	Catalog No.	Order Date	Order Status
Int Test Reg V173	RealTime ready Custom Panel 24*4	000000453	04743725001	2009-09-03	Evaluation
Syntenest US 002	RealTime ready Custom Panel 64*16	000000240	04720666001	2009-08-24	Production
Syntenest US 001	RealTime ready Custom Panel 24*4	000000231	04743725001	2009-08-24	

▼ My Ordered Single Assays

Assay ID	Gene Symbol (HGNC)	Alias	Configuration No.	Order Date	Order Status
Catalog Assay for ABCA1	ABCA1	ABC-1, ABC1, CERP, FLJ14958, HOLET1, MGC164864, MGC165011, TGO	000000259	2009-08-24	
Designer Assay for CCDC62	CCDC62	FLJ40344, TSP-NY	000000462	2009-09-03	Evaluation

Order Status

Evaluation

Production

Make a Copy

View Panel

Add to Shopping Cart

You can track the order status and retrieve the final content of your panels in form of a text file as soon as production of your panel is completed at Roche. This text file can be transferred to the sample editor of the LightCycler® 480 Software.


On the “My Orders” page, you can also re-order panels and single assays.

Note that any Designer Assays you have previously ordered will now be available as Catalog Assays, with updated catalog numbers. If you choose to re-order a previous Designer Assay without using the Configurator, for example, when ordering by phone, be sure to use the new catalog number.

To make ordering easier, we recommend using the RealTime ready Configurator to place your orders and re-orders. Using the Configurator ensures that the correct catalog number and price are automatically matched to your order.

Ordering Information

Product	Cat. No.	Pack Size
RealTime ready Custom Assays	Please go to www.realtimeready.roche.com	
RealTime ready Custom Panels	Please go to www.realtimeready.roche.com	
RealTime ready Focus Panels		
RealTime ready Human Reference Gene Panel, 96	05 339 545 001	2 plates (each containing 96 assays)
RealTime ready Human Reference Gene Panel, 384	05 467 675 001	2 plates (each containing 384 assays)
RealTime ready Human Apoptosis Panel, 96	05 392 063 001	2 plates (each containing 96 assays)
RealTime ready Human Apoptosis Panel, 384	05 339 316 001	2 plates (each containing 384 assays)
RealTime ready Human ABC Transporter Panel, 96	05 339 324 001	2 plates (each containing 96 assays)
RealTime ready Human ABC Transporter Panel, 384	05 467 713 001	2 plates (each containing 384 assays)
RealTime ready Human Cell Cycle Regulation Panel, 96	05 339 359 001	2 plates (each containing 96 assays)
RealTime ready Human Cell Cycle Regulation Panel, 384	05 467 683 001	2 plates (each containing 384 assays)
RealTime ready Human GPCR Panel, 96	05 353 068 001	2 plates (each containing 96 assays)
RealTime ready Human GPCR Panel, 384	05 467 705 001	2 plates (each containing 384 assays)
RealTime ready Human Nuclear Receptor Panel, 96	05 339 332 001	2 plates (each containing 96 assays)
RealTime ready Human Nuclear Receptor Panel, 384	05 467 691 001	2 plates (each containing 384 assays)
Real-Time PCR Reaction Mixes		
RealTime ready DNA Probes Master	05 502 381 001	5 x 1 ml (5x conc.)
LightCycler® 480 Probes Master	04 707 494 001	5 x 1 ml (2x conc.)
Additional Products Additional Products		
Transcriptor First Strand cDNA Synthesis Kit	04 379 012 001	1 Kit (50 reactions)
LightCycler® 480 Sealing Foil	04 729 757 001	50 foils



Disclaimer/trademark statement:

For life science research only. Not for use in diagnostic procedures.

REALTIME READY, LIGHTCYCLER and NIMBLEGEN are trademarks of Roche.

License Disclaimer:

A license to perform the patented 5' Nuclease Process for research is obtained by the purchase of (i) both Authorized 5' Nuclease Core Kit and Licensed Probe, (ii) a Licensed 5' Nuclease Kit, or (iii) license rights from Applied Biosystems.

This product is an Authorized 5' Nuclease Core Kit. Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155, 5,677,152, 5,773,258, 5,407,800, 5,322,770, 5,310,652, 5,210,015, 5,487,972, and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. Separate purchase of a Licensed Probe would convey rights under the applicable claims of US Patents Nos. 5,538,848, 5,723,591, 5,876,930, 6,030,787, 6,258,569, 5,804,375 (claims 1-12 only), and 6,214,979, and corresponding claims outside the United States. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

This LightCycler® 480 Real-Time PCR System is a real-time thermal cycler licensed for research use only under U.S. Patent No. 6,814,934 and corresponding claims in its non-U.S. counterparts, and under one or more of U.S. Patents Nos. 5,475,610, 5,602,756, 6,703,236, 7,238,517, 7,504,241, 7,537,377 or corresponding claims in their non-U.S. counterparts, owned by Applied Biosystems, LLC. No rights are conveyed expressly, by implication or estoppel to any other patent claims. The LightCycler® 480 Real-Time PCR System is not cleared for in vitro diagnostic use in the U.S. For further information on purchasing additional rights, contact the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California, 94404, USA.

No right is conveyed expressly, by implication, or by estoppel under any other patent claim, such as claims to apparatus, reagents, kits, or methods such as 5' nuclease methods. License rights to practice PCR methods under the foregoing Roche patents for research and other non-in vitro diagnostic applications may be purchased from Applied Biosystems or may be obtained by purchasing licensed reagents from Roche, Applied Biosystems, or other Authorized third party. For information on purchasing licenses for research and other non-in vitro diagnostic applications, contact the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

A license to practice PCR methods with real-time detection under patents of F. Hoffmann-La Roche Ltd, and Roche Molecular Systems ("Roche") for in vitro diagnostic applications may be purchased from Roche or may be obtained by purchasing in vitro diagnostic reagents from Roche or any other Authorized third party. The LightCycler® 480 Real-Time PCR System is not cleared for in vitro diagnostic use in the U.S.

Roche Applied Science offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our home page, **www.roche-applied-science.com**, and our Special Interest Sites including:

- The Universal ProbeLibrary:
<http://www.universalprobelibrary.com>
- RealTime ready Custom Assays and Panels
<http://www.realtimeready.roche.com>
- The LightCycler® System family for real-time, online PCR:
<http://www.lightcycler-online.com>
- The MagNA Pure System family for automated nucleic acid isolation:
<http://www.magnapure.com>
- Amplification – Innovative tools for PCR:
<http://www.roche-applied-science.com/pcr>
- DNA & RNA preparation – Versatile Tools for Nucleic Acid Purification:
<http://www.roche-applied-science.com/napure>

Published by

Roche Diagnostics GmbH
Roche Applied Science
Werk Penzberg
82372 Penzberg
Germany

© 2009 Roche Diagnostics GmbH.
All rights reserved.